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- Langdell RD, Wagner RH, Brinkhous KM. Effect of antihemophilic factor on one-stage clotting tests. *J Lab Clin Med* 1953;41:637-47.
- Turi DC, Peerschke EI. Sensitivity of Three Activated Partial Thromboplastin Time Reagents to Coagulation Factor Deficiencies. *Am J Clin Pathol*; 1986; 85: 43-49.
- Brandt JT, Triplett DA, Laboratory Monitoring of Heparin. Effect of Reagents and Instruments on the APTT, *AM.J.Clin. Path.* 1981; 76: 530-537.
- Van der Velde EA, and Poller L. The APTT Monitoring of Heparin - The ISTH/ICSH Collaborative Study. *Thromb Haemos*; 1995; 73: 73-81.
- Greaves M. Limitations of the laboratory monitoring of heparin therapy. Scientific and Standardization Committee Communications: on behalf of the Control of Anticoagulation Subcommittee of the Scientific and Standardization Committee of the International Society of Thrombosis and Haemostasis. *Thromb Haemost*; 2002 Jan;87(1):163-4.
- Kelsey PR, The Diagnosis of Lupus Anticoagulants by the Partial Thromboplastin Time. *Thromb. Haemost*; 1984; 52: 172-175.
- Tripodi A, Chantarangkul V, Asti D, et al. Activated partial thromboplastin time: results of a case-control study evaluating six commercial reagents in assessing the risk of venous thromboembolism. *J Thromb Haemost* 2006;4:1407-9.
- Tripodi A, Chantarangkul V, Martinelli I, et al. A shortened activated partial thromboplastin time is associated with the risk of venous thromboembolism. *Blood* 2004; 104,12: 3631-3634.
- Boekel E, Bartels Piet. Abnormally Short Activated Partial Thromboplastin Times Are Related to Elevated Plasma Levels of TAT, F1+2, D-Dimer and FVIII:C. *Pathophysiol Haemost Thromb* 2002;32:137-142.
- Hron G, Eichinger S, Weltermann A, Quehenberger P, Halbmayer WM, Kyrle PA. Prediction of recurrent venous thromboembolism by the activated partial thromboplastin time. *J Thromb Haemost* 2006;4:752-6.
- CLSI Document H21-A5, vol. 28, No. 5. Collection, transport, and processing of blood specimens for testing plasma-based coagulation assays and molecular hemostasis assays; approved guideline - 5th edition, 2008.
- CLSI Document H3-A6. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard - Sixth Edition, vol. 27, No. 26, 2007.
- Ray M. Stability of the activated partial thromboplastin time used to monitor unfractionated heparin. *J Thromb Haemost* 2008; 6: 1817-9.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. *Clin Chem* 1986;32:470-475.
- Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". *Eur J Clin Chem Clin Biochem* 1996;34:385-386.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. *Ann Clin Biochem* 2001;38:376-385.
- Stroobants AK, van Dam W, Bakker B, et al. Interference study of direct thrombin inhibitors and anti-Xa inhibitors on hemostasis assays on a cobas t 411 system, ISTH congress 2015, poster abstract.
- Cawkwell RD. Patient's age and activated partial thromboplastin time test. *Thromb Haemostasis* 1983;39:780-781.

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- 19 Gallistl S, Muntean W, Leschnik B, et al. Longer aPTT Values in Healthy Children than in Adults: No Single Cause. *Thrombosis Research* 1997;88:355–359.
- 20 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11):783-790.

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© 2015, Roche Diagnostics
Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim